

## Orris Breast Enhancement System

Effective and Proven Successful in Conditions like:

1. Sagging or shrinking of Breast with age.
2. Loosening of Breasts after pregnancy and feeding the child.
3. Weakness of mammary glands.
4. Very little or almost no development of Breasts due to lesser production of female hormones at puberty.
5. For the treatment of flat Breast area and flabby breast.
6. Depressed Nipples.
7. Female frigidity.
8. For beautifying yourself and becoming more attractive.
8. Lifting up and tightening the breast line.

## What is Orris Breast Enhancement System?

The Orris Breast Enhancement System is the only clinically proven, non-surgical alternative for breast enlargement and shaping, effectively and naturally. The Orris Enhancement System is safe, effective and painless. This system applies a gentle, sustained pull to the breast area resulting in the growth of new breast tissues. This process applies the principle of tension induced tissue growth. That has been widely used in medicine for over 30 years in various reconstructive procedures.

## The Orris Breast Enhancement System

This system is designed for women of all ages. Developed by medical specialists.

It consist of a powerful manually operated vacuum pump, two different sized domes made up of a special material called acrylic. These domes are fitted with specially engineered rubber rims for your comfort and ease of use



The dome is attached with the pump by a silicon tube and held in place after applying a lubricant. The negative pressure (Vacuum) is created by pump. This creates a gentle three-dimensional pull around the breast. This gentle pull stretches the cells, which respond by replicating or multiplying themselves. Until the resulting gap is filled. The out come of this exercise is new tissue growth and therefore larger and firmer breast giving you a new look.

## How Does This System Work?

The Orris Enhancement System works on natural principle that " which ever part of the human body is exercised results in GROWTH".


Pic. I

The baseline breast tissue consisting of sacks. These sensitive tissues having a natural tendency of expansion.




Pic. II

When activated by the Orris Vacuum System, the breast is gently pulled out. This stretches the sacked breast tissues.




**Pic. III**

*Due to the regular application of negative pressure the tissues start increasing their number to sustain the pull. They also become more firm.*



**Pic. IV**

*After regular exercise of a few days, the breast is enlarged and lifted by newly created healthy and more firm tissues. The regular exercise can add 2 to 3 cups or approximately 100cc to 300cc to the original breast size.*



**What are the side affects?**

*It is a medically proven and tested non-surgical device by FDA, USA which is based upon physio therapy principles. It is for external use and works on natural principle of " exercise and grow". Therefore it has no side affects.*

**What if my skin irritates?**

*The lubrication with gel, cold cream or any oil is recommended only to prevent skin irritation. In any case the control is in your own hands. Stop creating any more pressure if you start feeling even a slight pain.*

**Precautions**

*This device should not be used during pregnancy.*

*Don't transfer to others so as to avoid any possible skin disease.*

**Frequently Asked Questions**

**1) Who can use it?**

*Women of all ages and with all body types can use this system successfully and easily in the privacy of their own home.*

**2) How much time is required to see the results?**

*Every individual has her own level of body flexibility and growth rate. However, there can be a significant difference with in 3 to 6 weeks of regular use of five days a week.*



*Also Deals in men's impotency devices*